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
症状	マイルド	モデレート	ストロング
低緊張 回内 	Visible medial arch. Mild heel eversion and forefoot abduction. Can correct when prompted. Can be manually corrected with no resistance. PattiBob® 300-1225 in.	Reduced medial arch. Moderate heel eversion and forefoot abduction. Can improve when prompted. Can be manually corrected with mild resistance. Cricet® 400-900 in.	Absent medial arch. Strong heel eversion and forefoot abduction. Cannot improve when prompted. Can be manually corrected with moderate resistance. Chipmunk® 400-1225 in. JumpStart®/Leap Frog® 400-900 in.
高緊張 回内、回外 	HIGH TONE PRONATION Visible medial arch. Mild heel eversion and forefoot abduction. Can correct when prompted. Can be manually corrected with mild resistance. Chipmunk 400-1225 in. Pronation only	HIGH TONE SUPINATION Mildly increased medial arch. Mild heel inversion and forefoot adduction. Can correct when prompted. Can be manually corrected with mild resistance. JumpStart Leap Frog 400-900 in.	HIGH TONE PRONATION Reduced medial arch. Moderate heel eversion and forefoot abduction. Can improve when prompted. Can be manually corrected with moderate resistance. JumpStart Bunny 400-900 in. Wrap-around strap
遊脚期 異常歩行 	Lands heel first. No obvious compensations of the knee and hip. Occurs almost never (less than 80% of the time). Can control when prompted. JumpStart Bunny 400-900 in. Wrap-around strap	Lands foot-flat, accompanied by pronation/supination. Some compensations of the knee and hip. Occurs almost always (80% of the time). Can improve when prompted. JumpStart Bunny 400-900 in. Wrap-around strap	Lands forefoot-first, accompanied by pronation/supination. Marked compensations of the knee and hip. Occurs constantly (100% of the time). Cannot control when prompted. JumpStart Bunny 400-900 in. Wrap-around strap
つま先歩き 	Ankle plantarflexion: 0°. Occurs occasionally (less than 50% of the time). Can correct when prompted. Can be manually corrected with mild resistance. JumpStart Bunny 400-900 in. Posterior strap	Ankle plantarflexion: 0-2°. Occurs frequently (50% of the time). Can improve when prompted. Can be manually corrected with moderate resistance. JumpStart Bunny 400-900 in. Posterior strap	Ankle plantarflexion: 2° or more. Occurs constantly (100% of the time). Cannot correct when prompted. Can be manually corrected with strong resistance or cannot be corrected. JumpStart Kangaroo® 400-900 in.
膝の過伸展 	Gentle knee extension: 0-2°. Occurs occasionally (less than 50% of the time). Can correct when prompted. Can be manually corrected with mild resistance. JumpStart Bunny 400-900 in. Posterior strap	Marked knee extension: 2-5°. Occurs frequently (50% of the time). Can improve when prompted. Can be manually corrected with moderate resistance. JumpStart Bunny 400-900 in. Posterior strap	Significant knee extension: 5° or more. Occurs constantly (100% of the time). Cannot correct when prompted. Can be manually corrected with strong resistance. JumpStart Kangaroo® 400-900 in.
かがみ姿勢 	Gentle excess dorsiflexion and knee flexion: 5-10°. Occurs occasionally (less than 50% of the time). Can correct when prompted. Can be manually corrected with mild resistance. JumpStart Bunny 400-900 in. Wrap-around strap	Marked excess dorsiflexion and knee flexion: 10-15°. Occurs frequently (50% of the time). Can improve when prompted. Can be manually corrected with moderate resistance. JumpStart Bunny 400-900 in. Wrap-around strap	Significant excess dorsiflexion and knee flexion: 15° or more. Occurs constantly (100% of the time). Cannot correct when prompted. Can be manually corrected with strong resistance or cannot be corrected. JumpStart Bunny 400-900 in. Wrap-around strap
ポジショニング (姿勢管理) 	Accompanied by mild pronation/supination. Frequent assisted ambulation. Can be manually corrected with moderate resistance. JumpStart Bunny 400-900 in. Posterior strap	Accompanied by moderate pronation/supination. Occasional assisted ambulation. Can be manually improved. JumpStart Kangaroo 400-900 in.	Accompanied by strong pronation/supination. Assisted transfers only; or non-weight-bearing. Cannot be manually corrected. JumpStart Kangaroo 400-900 in.